

TEST SCOPE for all-inclusive ANALYSIS KITS

1	2	3	4	5
Standard Compressor	Standard Compressor	Standard Compressor	Standard Compressor	Standard Compressor
		OELCHECK recommendation Compressor Example lab-report		
<p>Test results</p> <p>Elements with ICP PQ index Viscosity 40 °C Viscosity 100 °C VI (viscosity index) Visual evaluation Oxidation Colour index IR index Water with FT-IR (%)</p>	<p>Test results</p> <p>Elements with ICP PQ index Viscosity 40 °C Viscosity 100 °C VI (viscosity index) Visual evaluation Oxidation IR index Water with FT-IR (%) Colour index Neutralisation Number (NN) or AN</p>	<p>Test results</p> <p>Elements with ICP PQ index Viscosity 40 °C Viscosity 100 °C VI (viscosity index) Visual evaluation Oxidation Neutralisation Number (NN) or AN Colour index Water acc. to Karl Fischer IR index</p>	<p>Test results</p> <p>Elements with ICP PQ index Viscosity 40 °C Viscosity 100 °C VI (viscosity index) Visual evaluation Oxidation Neutralisation Number (NN) or AN Colour index Water acc. to Karl Fischer Particle count ISO 4406, Cleanliness class IR index</p>	<p>Test results</p> <p>Elements with ICP PQ index Viscosity 40 °C Viscosity 100 °C VI (viscosity index) Visual evaluation Oxidation Colour index Water acc. to Karl Fischer Oxidation stability (RULER) IR index</p>

ANALYSIS KIT (all-inclusive)



- Original OELCHECK sample-bottle (100 ml) for analyses in the laboratory in Germany.
- Pre-addressed envelope.
- UPS return service (within Germany free of charge).
- Sample Information Form including an adhesive barcode label.
- Laboratory tests within 48 hours.

- Laboratory Report completed by a very detailed diagnostic statement of a mechanical engineer.
- Detailed representation can only be carried out at additional costs.
- Dispatch of the Laboratory Report via email, mail or fax.
- Password-safed online access to your laboratory reports and analysis data at www.lab.report.